**WP 3**

***D3.5 – User manuals***

**Category:** Deliverable  
**Author(s):** Heidi Similä, VTT  
**Date:** 25.5.2015  
**Status:** Final, v1.0

## Table of Contents

1	Introduction .....	3
2	Modules .....	4
2.1	Home exercising system .....	4
2.2	Activity monitor .....	5
	Appendix I: First version home system instructions for the end-user .....	6
	Appendix II: Second version home system instructions for the end-user .....	14
	Appendix III: Home system instructions for the Spanish testing - end-users .....	32
	Appendix IV: Installation guide for test coordinator – Spain .....	41

# 1 Introduction


This document gives a brief overview of user manuals created for different modules of the AiB solution. The user manuals include 1) usage instructions for the (older adult and professional) end-users and 2) installation guides for field trial organizers. Some of the manuals are written in the native language of the test users, e.g. in Finnish. Due to evolving nature of system development it was decided not to translate all the documents in English, since these manuals become obsolete with new versions of the system.

This document is targeted for project partners and professional users of the system.

## 2 Modules

### 2.1 Home exercising system

The home exercising system developed by VTT consists of a touch screen based exercising software and an accelerometer interlinked accelerometer software for automatic data download and analysis. The following picture illustrates the AiB development process. Different versions of the home software were delivered for the end-users, specifically version 1.0 at 10/2013 and version 2.0 at 05-06/2014.

Finnish instructions for the first version can be found in Appendix I ja for the second version in Appendix II. These instructions are targeted for the end-users, especially older adults using the system at home.

The home software developed by VTT was tested in Spain by health care professionals and older adults. For that purpose a end-user instructions were prepared in English (Appendix III) and also installation guide for local test coordinator (Appendix IV).

## 2.2 Activity monitor

The activity monitoring system may also be used separately apart from touch screen PC software. It consists of a USB accelerometer (x16-2 by Gulf Coast Data Concepts, LCC <http://www.gcdadataconcepts.com/x16-2.html>) and a USB docking station that is connected to the main home device (touch screen PC).

The related software components are:

- AiBactmon: Data download and analysis application
- Fileupload: Data upload client for data transfer to server

More detailed installation guide for data download and analysis software can be found in Appendix V.

# Appendix I: First version home system instructions for the end-user

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# Appendix III: Home system instructions for the Spanish testing - end-users

## AiB Spain setup user manual


V1.0

Heidi Similä

## Home screen

Note. button selection might be different in your setup.

# Home exercise

1. Push “Home exercise” button on the screen
2. Select type of exercise
3. Select level of exercise
4. Start exercising

This workout has 6 exercises

Stretching 3 **1/6**

Close window

Close

Stretching 3

1. Exercise

Exercise:

Sit on the front edge of a chair with the soles of your feet on the floor.

Straighten one leg in front of you with your heel on the ground. Bend forward. Feel the stretch in the back of the thigh muscles of the straight leg.

Keep stretching 5 seconds.  
Repeat 15 times.  
Repeat with the other leg.


Next

© ODL/PhysioTools Ltd

The screenshot shows a user interface for a home exercise routine. On the left, a woman is shown sitting on a chair, leaning forward with one leg extended. The text 'Stretching 3' and '1/6' are displayed, with '1/6' circled in red. An arrow points from the text 'This workout has 6 exercises' to the '1/6' indicator. On the right, there are two dropdown menus: 'Stretching 3' and '1. Exercise'. Below these, the exercise instructions are listed. At the bottom right, a 'Next' button with a right arrow is circled in red, with an arrow pointing to it from the text 'Go to next exercise'. At the top right, a 'Close' button with a close icon is circled in red, with an arrow pointing to it from the text 'Close window'.

Go to next exercise

Once all the exercises are done in the current workout, the system asks:

Save workout to the calendar

## Calendar: displays performed exercises

# Add own exercises to the calendar

E.g.:

Choose the date, sport, type and duration from the boxes, when you are ready tap save.

Date: 24. | October | 2013

Time: 9 | 0

Sport: Jogging

Type: Endurance

Duration: 0 | Hours | 50 | Minutes


Return

Add ✓

Save exercise to the calendar



## Adding questions/questionnaires to the system

- The questions appear as pop-up screens on **scheduled times**
- Questionnaires: If several questions are scheduled for the same time, they appear one after the other as a sequence

Sign in by touching the screen on the bottom right corner

Use password: a2e2admin →

1. Select "Events"

2. Select "Questions" for adding a new question

You can see the list of questions and occurrences currently inserted

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Question can be deleted here

Add new question here

1. Fill in the question and answer alternatives

New question

Question:  
How are you feeling today?

Answer 1:  
Adequate

Answer 2:  
Moderate

Answer 3:  
OK

Answer 4:  
Good


Answer 5:  
Excellent

Leave these blank

Occurrences | Ready

2. Add time when the question should appear


### List of times when the question appears

Add new occurrence when the question should appear


To see the questions: sign out from the admin system and the question should appear on the home view on scheduled times

E.g.:


If there is no internet connection, this screen appears

Note. The videos won't work without internet connection!

## Accelerometer

Can be used with the belt or in the pocket

# Accelerometer

- 1) Remove the accelerometer from the docking station
- 2) Place the “hat” on
- 3) Wear the sensor in the belt or in the pocket
- 4) Go for a walk, gym, running, golfing, ...
- 5) Place the sensor to the docking station when coming back
- 6) Step count is automatically uploaded to the calendar

## **Appendix IV: Installation guide for test coordinator – Spain**

- Removed from public Deliverable version